
Tosoh Bioscience

AIA® 360

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
MEASUREMENT OF SERUM PROGESTERONE IN MARES: A METHOD COMPARISON STUDY

C. Lubas, A. Gavazza, A. Rota, D. Panzani, S. Demi
Department of Veterinary Clinic - University of Pisa - Italy

INTRODUCTION
Serum progesterone (P4) concentration in the mare can be used to determine the phase of the estrous cycle and odds useful information in the optimal and maintenance of early pregnancy (1-4).
The aim of this study was to compare several methods/techniques to assess P4 values in mares, including fluorescence enzyme immuno assay (FEIA) (2 different assays), chemiluminescence (2 different assays), and liquid chromatography tandem mass spectrometry (LC/MS).

MATERIALS & METHODS
Forty serum samples from mares in different estrous cycle or pregnancy stages were collected and stored at -20°C until their testing. All samples were analyzed with the FEIA marketed by Tosoh Biosciences (TB), Belgium-Japan (AIA 360). Two types of investigations were arranged:
- Recovery study including low, medium and high values for equine P4 (1)
- Comparison study arranging set of 15 samples including low, medium and high values of P4 using the regression analysis (2). The methods/techniques chosen were: Chemiluminescence by IPC, Immulite (Siemens Medical); at two different commercial kits (Chem-1 and Chem-2); MEIA by ADVIA (Abbott Diagnostic); and ultra performance liquid chromatography tandem mass spectrometry (LC/MS) (Waters) (5).

The recovery study is reported in the table (6). In figures 1-4 are reported the comparisons carried out.

Horse P4	Rec. low (2.57 ng/mL)	Rec. medium (10.35 ng/mL)	Rec. high (20.61 ng/mL)	Mean recovery	Proportional error
Low value (2.71 ng/mL)	100%	102%	100%	100%	-26
Medium value (10.30 ng/mL)	100%	103%	101%	101%	-30
High value (19.88 ng/mL)	91%	95%	102%	96%	7

Recovery study for FEIA-TB showed results over 100% especially when testing low and medium equine P4 values. FEIA-TB had a very good agreement with the two chemiluminescence assays (2 > 95), while the comparison with a similar technique (MEIA-ADV) showed a lower (2 (3)). The lowest (2 (3)) was found between the FEIA-TB and LC/MS values. In all comparisons arranged, values obtained from FEIA-TB were looked slightly higher than those obtained from the other methods/techniques.

CONCLUSION
The P4 measurement using FEIA-TB showed good agreement when compared with different methods/techniques. It should be pointed out that the FEIA-TB gave results consistently slightly higher than other methods/techniques compared on non-specific reference range should be arranged for its clinical use in equine reproduction.

- Selected References
(1) Allen WR - Reprod Dom Anim. 36:121-131, 2001; (2) Elmore RG, Kleppe LH, Venzor DD, Meyers PJ - Vet Med. 83:294-297, 1988; (3) Nayak D, Henderson G, Belgodai J et al - Theriogenology. 61:200-214, 2004; (4) Squires EL, Westworth BC, Clifton DJ - J Am Anim Hosp. 19:276-287, 1979; (5) McDonald M, Mahony E, McElroy S - J Am Vet Assoc. 90:343-349, 2010; (6) <http://www.vetmed.purdue.edu/ans27.htm>, interference and recovery experiments - accessed Sept. 2011

Serum Cardiac Troponin Measurement by Fluorescence Enzyme Immuno Assay. A Comparison Study and its Application in Canine Practice

A. Gavazza, L. Carlucci, S. Demi, G. Lubas
Dept. of Veterinary Sciences, University of Pisa, Italy


Introduction
Measurement of serum Cardiac Troponin (cTn) concentration in dogs to assess myocardial damage or injury and as prognostic indicator of future cardiovascular disease has been used. Currently, various cTn assays are available producing different results with, consequently, using different critical cut-off values. Therefore, cTn values obtained from different assays cannot be interchanged for cut critical values.

Materials and Methods
Forty-seven serum samples collected from healthy dogs (assessed by physical exam and complete blood count, biochemistry, and urinalysis) were used to set the Reference Range (RR) of cTnI FEIA-TB (Fluorescence Enzyme Immuno Assay Tosoh Biosciences). These samples were also tested with cTnI Chem (Chemiluminescence Immunoassay). In addition, 147 serum samples collected from MVD (Mitral Valve Disease)-affected dogs and staged according to Cardiovascular Disease, ACVIM Consensus Statement (stage A not studied, B #61, C #10, and D #51) were tested both with cTnI FEIA-TB and cTnI Chem. Statistical analysis was arranged for data collected (RR, correlation, r and regression β , Spearman, rho for dogs in different CVD-stages).

Results
RR for cTnI FEIA-TB were set at 0.0-0.09 ng/mL. In comparison, RR for cTnI Chem were 0.05-0.24 ng/mL.

dogs	FEIA-TB > RR	Chem > RR	r	β	rho
healthy	0	0	-	-	-
B	61	12	0.99	0.99	0.72
C	19	14	0.96	0.96	0.97
D	11	9	0.95	0.90	0.93

Comparison of the RR for cTnI FEIA-TB assays (dogs) tested by reference ranges

Conclusion
RR of cTnI with FEIA-TB in healthy dogs are low in comparison with those reported for cTnI Chem but they match very closely to identify dogs with myocardial damage. Indeed, a cut-off value in good agreement when the comparison is carried out in CVD-stages of MVD-dogs as FEIA-TB correctly identify myocardial damage.

References are available at: agavazza@vet.unipi.it

ACCURATE DETERMINATION OF SERUM PROGESTERONE USING A FLUORESCENCE ENZYME IMMUNOASSAY IN THE BITCH

A. Monino¹, P. Dri¹, M. Beccaglia², G. Milite³

- 1- Clinica "Animal Care", Via Maù 18, 33035 Marignacco (UD), Italy;
2- Ambulatorio Veterinario Beccaglia, Via Volta 2, 20851 Lissone (MB), Italy;
3- Scientific Consultant, Via Susans 40, 33100 Udine, Italy. g.iuspaula.milite@justvetnet.it

International Symposium on Canine and Feline Reproduction

ISCFR 2012 - EVSSAR 2012

26 - 29 July 2012
WHISTLER, B.C. CANADA

A comparison of two different Enzyme-Linked Fluorescent Assay Instruments for canine progesterone assay in the clinical practice

Majolino G, Alonge S, Benin D, Conti C
1) Ambulatorio Veterinario Dr. R. Ranieri - Dr. G. Majolino, Collecchio, Parma, Italy
2) Società Veterinaria "Il Molleggiato" Srl, Sesto Candelo, Varese, Italy
3) Scientific Consultant Services, Ferrara, Ferrara, Padova, Italy

Introduction
Progesterone (P4) concentration in peripheral blood is considered the best method to detect the best fertile period, to predict the whelping day, to control luteal insufficiency, to verify luteolysis prior to parturition in the bitch (1). Mini Vidas (Biomérieux, Marcy-l'Étoile, France) system allows to obtain rapid and reliable results for determination of P4 by the enzyme-linked fluorescent assay (ELFA) with high agreement (mean deviation 15%, correlation coefficient 0.989) if compared to a radioimmunoassay (RIA) validated in the dog (2).

Moreover, in our clinics we processed 1500 progesterone analyses yearly by MiniVidas: 90% of bitches monitored during estrus whelped 63 days (+/-) after the estimated ovulation day, as in the literature (1).
Aim
To determine the correlation between two different Enzyme Immuno Assay analyzers for P4 measurement: Tosoh AIA360 - fluorescent enzyme immunoassay method (FEIA) (3) and MiniVidas - ELFA method (4) and to calibrate the FEIA system to the ELFA one, in order to obtain results interchangeable in everyday clinical practice.

Both groups included different breeds, size and age bitches:
1) for the detection of ovulation (58/106-52.46% and 47/103-44.47%)
2) with normal pregnancy (5/106-4.7% and 11/103-10.66%)
3) pregnant with hypoluteinism (3/106-2.84% and 5/103-4.85%)

Results
Group A (after regression between techniques (r=0.967), but mean value difference between results: [4.3625.52 (mean±SD)])
Group B (after regression analysis and calibration) High correlation (r=0.98) and lower mean value difference between results: [1.11.27 (mean±SD)]

Conclusions
Very high accordance level between the 2 different Enzyme Immuno Assay analyzers for P4 measurement after calibration. Measurements of progesterone for the detection of ovulation as well as in pregnant bitches with normal pregnancy or suffering from hypoluteinism obtained by both systems provided meaningful results. Correlating results on serum P4 concentration obtained by the 2 techniques indicate that they are interchangeable leading to the same clinical decision.

References
(1) Alivonoti et al., 2011, Breeding Management and Artificial Insemination of the Bitch in Canine and Feline Reproduction, 9th Edition, Pp 40-65
(2) Alonge et al., 2012, Quantitative Determination of progesterone (P4) in canine blood serum using an enzyme-linked fluorescence assay. Reprod Domest Anim, 37:40-47
(3) Lubas et al., 2014, Laboratory and Clinical Evaluation of a FEIA Method for Canine Serum Progesterone Assay. Reproduction in Domestic Animals, 39:18-24
21st EVSSAR CONGRESS 2012
EVSSAR 2012
June, 21st-23rd 2012
Vancouver, Canada